

Seat2Go

A lightweight soft padded positioning seat for children that need support while sitting on the floor, at the table or at play. It is simple to use in a home, school or when travelling.



drive

Your Seat2Go can be used as a floor sitter, booster chair and an insert in a stroller. They are suitable for some children that need support while sitting on the floor at the table or at play. The CRS2000 has a weight limit of 50lbs (23kgs) and the CRS3000 has a limit of 75lbs (34kg)



ADJUSTMENTS

The straps below (A-D) need to be adjusted to suit each user. The straps should be taut and supportive, but ensure that they are not too tight to restrict the user's breathing or comfort.


TO USE AS A FLOOR SITTER

1. Place the Seat2Go on a flat surface.
2. Open the locking clip (A) on the pelvic belt. Detach the two shoulder straps (B) from the velcro straps on the Seat2Go back. Open the locking clip (C) on the seat support (D) and pull the \ support apart.
3. Place the child in the Seat2Go and fasten the pelvic belt (A).
4. Wrap the chest support (D) snugly around the child. and secure using the velcro on the chest support and the locking clip (C).
5. Determine the correct fit for the shoulder straps (B) and attach the Velcro strips on the Seat2Go back. The shoulder straps should be supportive and taut but not pulling up on the chest support or uncomfortable for the child.



TO USE AS A BOOSTER SEAT

1. Remove the fastening belts from the pocket on the seat bottom (E) and Seat2Go back (F).
2. Place the Seat2Go on a suitable chair or in a stroller, wrap the belts around for the back (G) and seat (H) of the chair/stroller and secure by closing the locking clips.
3. Ensure that the chair is secured safely in place and that the belts are taut.
4. Insert the child as per the instructions for the floor sitter.



Custom Technologies Limited
Taupo, New Zealand

Tel: +64 7 378 5335
Fax: +64 7 378 5408

E: ctl@xtra.co.nz
W: customtech.org.nz
W: customtech.kiwi