

Adjusted by:

E-mail:

Adjusted for:

Tel:

Date:

Total weight (g):

The weight changes based on the choices made below.

U3-light has a rigid rear axle. To find a balance that suits you, there are two versions, **U3 light** or **U3 light L**. L is more easy balanced and easier to drive on the rear wheels. You can tune the balance by placing the backrest in three different position. U3 Light has a seat height of 47 cm in front and 43 cm at the rear without cushion.
Max userweight is 100 kg for seatwidth 33-45.

Save the guide to the "desktop" and open it with **Adobe Reader**, for best function Design your wheelchair by clicking in the circles in this guide. **Blue, underlined text is clickable to show info** or image.
BB = seatwidth, **HH**= backrest height

Original design

Alternative options

Other options

[info](#)

1) Chassis colour

- White** - Pearl White
- Black** - Anodic Black
- Clear

Optional colour - A higher price and longer delivery time for optional colour

[Optional colour](#)



2) Seat width (BB)

[article number](#)

- 33 cm
- 33 cm L
- Clear
- 36 cm
- 36 cm L
- 39 cm
- 39 cm L
- 42 cm
- 42 cm L
- 45 cm
- 45 cm L

Art.no includes: chassis, seat, calf strap and a seat bag



3) Backrest

Height (HH)

- U3 light** 523BBHH L
- Clear
- 25 cm
- 30 cm
- 35 cm
- Clear

Backrest upholstery Comfort

- U3 Light backrest with S3 upholstery - more padding and stronger ribbons than std - slightly heavier



4) Balance center of gravity (COG) [info](#)

- Position A** Least easy balanced
- Position B** Medium balance
- Position C** Most easily balanced
- Clear

The positioning of the backrest on the chassis affects the balance / COG of the wheelchair. "easy balanced" = easy to practice rear wheel balance. This can later easily be changed.



5) Castor

- 90 mm X** 2100222 (x2 / wheelchair)
- Clear

- 90 mm S3** 2100030 x2
- 120 mm S3** 2100020 x2



6) Footrest

- Footrest X** 36004BB
- Clear

- Footrest, Titanium** 36915BB
- Footrest extended, reinforced** 3606BB [needed for Freewheel](#)
- Plastic plate for foot bow** 36010BB

[Other footrest](#)

[List of options art.](#) _____



7) Brake

- One hand brake, right** [article number](#)
- One hand brake, left**
- Clear

- Brake S3** 4630000



8) Rear wheel - without pushrim

- 24" Spinergy X white spokes** 1100184 (x2 / wheelchair)
- 24" Spinergy X black spokes** 1100141 (x2 / wheelchair)
- Clear

- 25" Spinergy X white spokes** 1100185 x2
 - 25" Spinergy X black spokes** 1100146 x2
 - 26" Spinergy Spox black spokes** 1100100 x2
- [Rear wheel sizes](#)

[Other rear wheel](#)

[List of options art.](#) _____



9) Pushrim

- Aluminum X** 1904404 (5) (x2 / rullstol) 24" (25")
- Clear

- Titanium** 190DD11 x2
 - Paragrip** 19051DDBTx2
 - Tetragrip** 19052DDBTx2
- DD** = Wheel diameter

Maxgrip 24" (contains latex) 1903060 x2



10) Tire

- Schwalbe One 24", 25"** standard [art.no.](#)
- Clear

- Schwalbe **Right Run** 24",25",26"
- Rugged 24"**

Schwalbe **Marathon Plus** 24",25",26"

Accessories



Side guard

- [U2 light](#) [info](#)
article number
- [Side guard S3](#)
article number
- Clear



Side guard

- [Carbon, small cover](#)
article number
- [Carbon, large cover](#)
article number



Arm support

- [Arm support S3](#) [Video](#)
article number

Makes the wheelchair 15mm wider



Cushion

- [2,5 cm](#)
572BB40
- [5 cm](#)
575BB40
- Clear



Push handles

- [Foldable S3](#) [Video](#)
4245000
- Clear



Other

- [Freewheel](#) [info](#)
2600102
- [Tube protection](#)
9900251 (x3 / wheelchair) ___ pcs
- [Back rest wedge](#)
4710010 ___ pcs

Set-up adjustments



Backrest angle

Note measure (A) for desired back rest angle.
10-20 mm. Standard=18 mm

A= _____ mm



Footrest height

Note measure (A) for desired footrest height.
Standard = 70 mm

A= _____ mm

Extra