

Adjusted by:

E-mail:

Adjusted for:

Tel:

Date:



Panthera X has the chassis and backrest made in carbon fibre and a rigid rear axle.
The balance is adjusted by positioning the backrest in **4 different positions**.
Panthera X is delivered without anti-tip device and is intended for advanced users only.
Max user weight is **100 kg for seat width 33-45 cm**.

Save the guide to the "desktop" and open it with **Adobe Reader**, for best function
Design your wheelchair by clicking in the circles in this guide.
Blue, underlined text is clickable to **show info** or image.
BB = seatwidth, HH= backrest height

Original design

Alternative options

Other options

[info](#)



1) Seat width (BB) Chassis [art.no.](#) Q350BB11G [BB= seat width]

- 33 cm
- 36 cm
- 39 cm
- 42 cm
- 45 cm
- Clear



2) Backrest height 22-35 cm (HH) Q342BBHH [BB= seat width HH= backrest height]

- 22 cm
- 23 cm
- 24 cm
- 25 cm
- 26 cm
- 27 cm
- 28 cm
- 29 cm
- 30 cm
- 31 cm
- 32 cm
- 33 cm
- 34 cm
- 35 cm
- Clear



3) Balance [info](#)

- Position 1**
Most easy balanced
- Position 2**
Less easy balanced
- Position 3**
Medium balanced
- Position 4**
Least easy balanced
- Clear

The positioning of the backrest on the chassis affects the balance / COG of the wheelchair. "easy balanced" = easy to practice rear wheel balance. This can later easily be changed.



4) Castor

- 90 mm X**
Q2305130 (x2 / wheelchair)
- Clear

- 90 mm S3**
Q2100201 x2
- 90 mm for 26" rear wheels**
Q2305132 x2

- 120 mm for 24" rear wheels**
Q2305133 x2
- 4 x 1.5 in Soft roll caster x2**
Fork Q230513(4) (5) (6)
Ultra Caster 9002846 x2



5) Footrest

- Footrest Alu X**
Q36004 BB
- Clear

- Footrest Titanium**
Q36915BB
- Footrest Titanium, extended, reinforced**
Q36006BB [needed for Freewheel](#)
- Plastic plate for foot bow**
Q36010BB



7) Brake [art.no.](#) Q462BB02

- One hand brake, right**
- Clear

- One hand brake, left** Q462BB01



8) Rear wheel - without pushrim

- 24" Spinergy X black spokes** Q1100141
- 24" Spinergy X white spokes** Q1100184
- Clear [art.no.](#)

- 25" Spinergy X black spokes** Q1100146
- 25" Spinergy X white spokes** Q1100185
- 26" Spinergy X black spokes** Q1100176
- 26" Spinergy X white spokes** Q1100186

Art.no are for wheels with Schwalbe One

[Rear wheel sizes](#)



9) Tire

- Schwalbe One 24"**
- Clear

- Right Run 24",25",26"**
- Marathon Plus 24",25",26"**

- Rugged 24"**

[info](#)



10) Push rim

- Aluminum X**
Q1904404 (5) (x2 / wheelchair) 24" (25")
- Clear

- Titanium** **DD = Wheel diameter in inches**
Q190DD11 x2
- Paragrip**
Q19051DDBTx2
- Tetragrip**
Q19052DDBTx2

- Maxgrip 24" (contains latex)**
Q1903060 x2

Accessories




Side guard

- [U2 light](#)
arti.no.
24-25" wheels Q4420840
26" wheels Q4420843
- [info](#)
- Clear



Side guard

- [Carbon, small cover](#)
arti.no.
24" Q4420740
25-26" Q4420743



Side guard

- [Carbon, large cover](#)
arti.no.
24-25" Q4420840
26" Q4420743



Cushion

- [X 5 cm](#) seat depth 40 cm
Q47501BB
- Clear



Arm hook

- [Arm hook X](#)
Q4210100
- Clear



Transport protection

- [Housing](#)
Q9024009
- [Rear end](#)
Q9024012
- [Rear axle](#)
Q9024008



Other

- [Freewheel](#) [info](#)
Q2600102
- [Tube protection](#) _____ p
Q9900252 (x3 / wheelchair)

Folding and Locking Side guards

- Q549269

Set-up adjustments



Backrest angle

Note measure (A) for desired back rest angle
10-20 mm. Standard=18 mm

A= _____ mm



Footrest height

Note measure (A) for desired footrest height
Standard = 70 mm

A= _____ mm

Extra