



# Emys

## Technical sheet

English





## Main features

### Comfort

Each pad is made of foam, with a removable and washable covering realized with technical material for maximal comfort.

### Design

Minimal and soft lines highlight an hi-tech design blending beauty and ergonomics.

**New:** discover the high-visibility stripes for your safety.



### Fresh feel


The space between the slats and the specific materials used for the pads, make the backrest highly breathable.

### Flexibility

The central armonic spring steels and the modular system, allow the backrest to flex giving freedom of movement to the user.

**New:** with EMYS you can decide the level of flexibility of the backrest!


## Height






| 3V | 4V | 5V

- | 3 Vertebrae 27 cm
- | 4 Vertebrae 35 cm
- | 5 Vertebrae 43 cm

## Size



| M | L | XL |

-  | **M** for wheelchairs from 35 to 40 cm
-  | **L** for wheelchairs from 38 to 43 cm
-  | **XL** for wheelchairs from 40 to 45 cm

## Depth



**STANDARD**

5,5 cm



**CONTOUR**

9 cm



**DEEP-CONTOUR**


14,5 cm

Aluminum is the most adaptable material to make each slide of the backrest suit the user's body!










# Adjustable clamp for rigid frame wheelchairs





## Features:

- 1 height adjustment
- 2 rotation adjustment
- 3 depth adjustment
- 4 misure 19-22cm
- 5 easy unhook
- 6 safety rubber band











# Adjustable clamp for folding wheelchairs



### Features:

- 1 rigidizer bar to specify:  
**M:** 35 - 41 cm    **L:** 39 - 45 cm
- 2 hose clamp: 19-21 / 22-25 cm
- 3 height adjustment
- 4 depth adjustment
- 5 rotation adjustment
- 6 easy unhook
- 7 safety rubber band






## Finishing



**"A" bent slat**  
for a **more comfortable support**.

### Upper pocket pads

for a greater adherence to the slats and to **intensify the softness** of the backrest.



### Balance points

on the rigidizer bar to simplify the balance of the clamp for folding wheelchairs.

**2 rubber tubes** (internal diameter 22mm)  
easy adaptables to every wheelchair **to protect the backrest**.




### Highly breathable pads

removable and **washable** at 30°.



Choose your **style.**





Questions?  
Info?  
Help?

Visit our website:  
[www.tartadesign.it](http://www.tartadesign.it)

or contact us at:  
[info@tartadesign.it](mailto:info@tartadesign.it)



## **Tarta Design srl**

Via Cottonificio, 56  
33037 Pasiàn di Prato  
Udine - Italy  
T +39 0432 722 533  
F +39 0432 542 920

[www.tartadesign.it](http://www.tartadesign.it)  
[info@tartadesign.it](mailto:info@tartadesign.it)